Differential Glycoprotein Profiles in Severe and Asymptomatic COVID-19

Chad Pickering¹; Prasanna Ramachandran¹; Gege Xu¹; Rachel Rice¹; Srinivasa T. Reddy²; Faheem Guirgis³; Taia T. Wang⁴; Daniel Serie¹; Xin Cong¹ ¹InterVenn Biosciences, South San Francisco, CA; ²University of California, Los Angeles, CA; ³University of Florida, Gainesville, FL; ⁴Stanford University, Stanford, CA

Introduction

Glycosylation is the most common post-translational modification of proteir critically affecting their structure and function. We analyzed abundant serur and plasma proteins in 50 individuals with severe COVID-19 and in 22 subjec who had an asymptomatic course of COVID-19 using liquid chromatography mass spectrometry (LC-MS) for high-resolution site-specific quantification of glycopeptides coupled with high-throughput artificial intelligence/neural network (AI)-powered data processing of protein glyco-isoform distributions. As additional comparisons, we included 12 individuals seropositive for commoncold coronavirus (CC), 16 patients with bacterial sepsis, and 15 healthy subjects without history of coronavirus exposure. Data on 597 peptides and glycopeptides were analyzed using supervised and unsupervised machine learning (ML) techniques.

Methods

Peptides and glycopeptides from blood serum and plasma were identified in a prior study with data-dependent acquisition. Among them, a panel of >500 glycopeptides was selected for the targeted quantitative analysis of trypsindigested serum and plasma samples using the workflow below.



ns,	
m	
cts	
/ and	
of	



15

Color Key -4 -2 0 2 4 Row Z-Score



The ratio of **glycopeptide**

abundance relative to its

suitability glycopeptide

conjugate non-glycosylated

peptide is assessed in pooled

serum replicates by run order.

These plots feature five system

biomarkers in each abundance



LC/MS Quality Control Metrics

11



Healthy controls



Results and Discussion



34 biomarkers show statistically significant differences separating symptomatic COVID-19 from all other groups, FDR<0.05 Symp • Healthy • CC



Features extracted at the glycosylation and protein levels clearly differentiate symptomatic COVID-19 from asymptomatic/mild COVID-19 and other comparison cohorts such as bacterial sepsis and common cold (CC). A number of glycoproteomic biomarkers were shared between patients suffering from symptomatic COVID-19 and bacterial sepsis, likely representative of a systemic inflammatory response.

A subset of 34 glycoproteomic biomarkers were specific to severe COVID-19, potentially providing insights into specific pathological consequences of, or a particular predisposition for, a severe course of the infection.



	1	2	3
COVID+ symp	0	3	47
Other	49	0	0
Sepsis	0	14	2

Conclusions

Glycoproteomic profiles in patients with a severe course of COVID-19 infection are statistically significantly different from those who experienced an asymptomatic course. InterVenn's proprietary perspectIVTM platform, based on a high-throughput LC-MS/AI-ML approach, interrogates the blood glycoproteome at an unprecedented depth and scale, revealing a whole new domain of biology for the identification of novel biomarkers in a broad range of applications.